

Mercury analyzer
mercur



mercur

- mercur**
- PC controlled fully automatic Hg-analyzer, based on cold vapour technique and atomic fluorescence
 - Measuring range: low ng/L - mg/L-range

mercur plus

- PC controlled fully automatic Hg-Analyzer
- Based on cold-vapour technique and atomic fluorescence
- 2 gold collectors for single or double enrichment
- Measuring range: 0,05 ng/L to mg/L- range:

Measuring principle

- Cold vapour atomic fluorescence-spectrometry with excellent sensitivity, selectivity and linearity for Hg-analysis
- Hg vapour generation based on cold vapour technique (CV) with SnCl₂ as reducing agent
- Detection of the mercury fluorescence radiation at 253,7 nm at 90° angle to the direction of radiation source

Optical System

Spectrometer Atomic fluorescence spectrometer with focussing of the fluorescence radiation to the detector

Light source Hg-low pressure lamp with high intensity

Fluorescence cell Double sided metallized quartz flow cell, Dimensions 10 mm x 10 mm x 32

Detector Photomultiplier

System

Operation Mode Time-controlled flow injection, without or with autosampler, amalgamation unit as option

Gas-Liquid Separator Optimized gas-liquid separator, minimized foam formation and rinsing time

Gas-dryer Perma-Pur membrane dryer with optimized efficiency through counter flow principle

Gas- and Liquid transport

- 4-channel peristaltic pump with separated transport of reducing agent, acid and sample waste
- Separate low noise 1- channel peristaltic pump for sample only, saves sample volume and reagents

Reagents

- Reducing agent: 2% - 10% SnCl₂
- Acid: HCl 2,0% (Hg free)
- Reagent consumption: 2 mL Reducing agent and acid per measurement

Gas

- Argon (inert), Consumption max. 1L/Minute
- Input pressure 500kPa (5bar)

Analytical Parameter	
Detection Limit	$\leq 1 \text{ ng/L}$ (without enrichment) $\leq 0,1 \text{ ng/L}$ (with enrichment)
Dynamic measuring range	5 concentration decades ($1 \text{ ng/L} - 100 \mu\text{g/L}$) without electronic sensitivity control Carry over $< 0,5\%$ relative
Sample consumption	Approximately 1mL/measurement
Analysis time	Approximately 40 sec without amalgamation, approximately 100sec with amalgamation
Sample types	Liquid samples
Accessories	
Autosampler	AS-F / AS-FD – Autosampler with glass vials for fully automatic and contamination-free analysis
Other technical data	
Software	WinAAS®
Dimensions (W x H x D)	600 mm x 350 mm x 490 mm (without PC)
Weight	37 kg (without PC)
Ambient conditions	<ul style="list-style-type: none"> ▪ Corrosion-resistant to the analyzed samples ▪ Humidity max. 90%, non-condensing ▪ Temperature $+10^{\circ}\text{C}$ up to 35°C
Power Supply	230V ($\pm 10\%$); 50/60Hz, Fuse 3,15A, max. 600 VA power consumption
Technical standards	<ul style="list-style-type: none"> ▪ EN 61010-1-1; EN 61010-2-061; IEC 61010-2-061; EN 50082; EN 55011; EN 61326 ▪ ISO 9001; Certification EMC, CE Mark



Subject to changes in design and scope of delivery as well as further technical development!